

SCIETECHQMS

Computer System Validation Report

Installation · Operational · Performance Qualification
(IQ / OQ / PQ) · GAMP 5 · 21 CFR Part 11

SAMPLE DOCUMENT — FOR DEMONSTRATION ONLY



DOCUMENT TITLE

Computer System Validation Report

SYSTEM

ScieTechQMS v1.1.0 — Edge Intelligence Edition

ENVIRONMENT ID

ENV-SAMPLE-9F2A1C8B4E7D0031

DOCUMENT ID

STQMS-CSV-A3F1C92E7B4D

CHAIN OF TRUST (DB)

a3f1c92e7b4d8e0156c0aa91f2e3b4c5d6e7f8091a2b3c4d5e6f708192a3b4c5

REPORT DATE

2026-08-02

CLASSIFICATION

Confidential — Controlled Document (Sample)

OVERALL STATUS

PASSED

Table of Contents

1. Executive Summary	3
2. Installation Qualification (IQ)	4
3. Operational Qualification (OQ)	5
4. Performance Qualification (PQ)	6
5. Signature Log (21 CFR Part 11)	7

1. Executive Summary

OVERALL VALIDATION STATUS: PASSED

This sample Computer System Validation Report illustrates the IQ/OQ/PQ package produced by the ScieTechQMS CSV Engine for Edge Intelligence Edition v1.1.0.

Verification is performed by the on-device Chain of Trust engine. The database integrity hash (SHA-256) is embedded in PDF metadata to support tamper detection of the report package.

Automated GAMP 5 Validation Reports generated by our internal CSV Engine — local-first, no regulated data leaves the workstation.

Control	Result
21 CFR Part 11	COMPLIANT
ALCOA+	COMPLIANT
Audit Trail	INTACT
Chain of Trust Hash	a3f1c92e7b4d8e0156c0aa91f2e3b4c5d6e7f809...
Lockdown	CLEAR

2. Installation Qualification (IQ)

ID	Test	Evidence / Observed	Result
IQ-01	Operating System / Platform	Windows / Linux workstation (sample)	PASS
IQ-02	Host memory (RAM)	16 GB (host reported)	PASS
IQ-03	Disk / application data store	Local SQLite + app data accessible	PASS
IQ-04	Rust / Tauri runtime	Embedded backend responding (Tauri 2)	PASS
IQ-05	Local AI Engine (Phi-3 Mini)	Ready · SHA-256 8a83c7fb9049a9b2e922...	PASS
IQ-06	Environment ID (Machine Hash)	ENV-SAMPLE-9F2A1C8B4E7D0031	PASS

3. Operational Qualification (OQ)

ID	Test	Evidence / Observed	Result
OQ-01	Audit Trail Integrity Scan	128 entries · no anomalies (sample)	PASS
OQ-02	Chain of Trust verification	Hash chain valid · 128 blocks	PASS
OQ-03	Electronic Signature functionality	12 signatures · Part 11 controls active	PASS
OQ-04	Security Encryption (AES-256)	Local datastore cryptography verified	PASS
OQ-05	ALCOA+ / GxP compliance flags	COMPLIANT	PASS

4. Performance Qualification (PQ)

ID	Test	Evidence / Observed	Result
PQ-01	Database response time	Mean 12 ms · Max 14 ms (4 samples)	PASS
PQ-02	Data retrieval consistency	Consistent Chain of Trust hash across reads	PASS
PQ-03	Uptime / operational window	1d 0h 0m	PASS
PQ-04	Audit stream continuity	128 audit entries available for retrieval	PASS

5. Signature Log

Electronic approvals per 21 CFR Part 11. Sample signatures shown for demonstration.

Role	Name / Title	Date	Meaning of Signature
Prepared by	CSV Engineer — Sample	2026-08-02	I prepared this validation report and attest that the recorded evidence was generated by the ScieTechQMS Chain of Trust engine.
Reviewed by	QA Reviewer — Sample	2026-08-02	I have reviewed this validation report for completeness and technical accuracy against the approved CSV protocol.
Approved by	QA Manager — Sample	2026-08-02	I approve this Computer System Validation Report for release. The system is fit for intended GxP use.